

Quantifying an unknown concentration by calibrated titration
and integration of NMR peak areas
(YB, June 2003)

(1) The NMR tube contains an unknown amount of substance U to be quantified. The volume of material in the NMR tube will need to be precisely known (V). The spectral peaks of U must already be assigned and must be well resolved from any other peaks in the NMR spectrum. The NMR signal must be acquired under quantitative conditions.

$V_o \rightarrow$ known initial volume

$U_o \rightarrow$ unknown amount to be quantified

(2) A separate solution of the material U (or a suitable non-interacting reference) of known concentration ($C = U/V$) should be prepared.

(3) Acquisition hints for quantitative integrations:

Measure the longest T_1 and use $5 \cdot T_1$ for the recycle delay

Use a large sw to get a flat baseline

Center your data in the middle of the spectrum (o1)

Acquire a large number of dots (TD)

Make sure the signal-to-noise ratio is high (NS)

Use an optimal AQ time (fid going to zero at $\sim 1/3$ of the screen)

For decoupled measurements of X nuclei – you MUST avoid any NOE effects

Repeat for a few different independent measurements

(4) THE EXPERIMENT

i. Acquire an FID, FT, and integrate the original sample containing U_o/V_o . If you expect that after titration the signal will be truncated, then reduce the receiver gain at this stage, BEFORE acquiring the reference FID.

ii. add into the NMR tube an aliquot of known volume (V) of known concentration (C) so you now have $(U+U_o)/(V_o+V)$. Acquire an FID, FT, and integrate – using all the same parameters of step 4*i* (*i.e.* do NOT do an rga . . .)

iii. repeat step 4*ii* for several aliquots of material. A minimum of three aliquots should be used to make the calibration curve.

(5) Processing hints for quantitative integrations:

Do not use a window function, just do a straight fourier transform

Be sure that the baseline is flat on both sides of the peaks to be integrated

Include the ^{13}C satellites in the integration

Adjust the slope and bias of the integration

Use the same limits of integration for all spectra (wmisc \rightarrow intrng, rmisc \rightarrow intrng)

Use the same scale for all the integrations (“lastscal”)

Use deconvolution software, if lines are not resolved

(6) To create a calibration curve, start with a table summarizing your NMR results:

expno	U (μg) added	V (μl) added	V_o+V (μl)	conc of added (C')	Normalized peak area	Peak area of added
1	none	none	500.0	none	1.00	none
2	100	1.5	501.5	0.2 ($\mu\text{g}/\mu\text{l}$)	1.50	0.50
3	200	3.0	503.0	0.4 ($\mu\text{g}/\mu\text{l}$)	3.88	2.88
4	300	4.5	504.5	0.6 ($\mu\text{g}/\mu\text{l}$)	6.02	5.02

Then plot the peak area due to the added material versus the concentration added – this is your calibration curve.

Using the peak area of the material at its initial concentration equal to 1.00, the corresponding concentration can be determined directly from the graph.